
FISSA Documentation

Release 0.6.1

**Sander Keemink
Scott Lowe**

Mar 12, 2019

Contents:

1	FISSA	1
1.1	Usage	1
1.2	Installation	2
1.3	Folder Structure	3
1.4	Citing FISSA	4
1.5	License	4
2	API	5
2.1	fissa package	5
3	Contributing code	17
3.1	Reporting Issues	17
3.2	Documentation	18
3.3	How to contribute	18
3.4	Notes	20
4	Changelog	21
4.1	Version 0.6.1	21
4.2	Version 0.6.0	21
4.3	Version 0.5.3	22
4.4	Version 0.5.2	22
4.5	Version 0.5.1	22
4.6	Version 0.5.0	22
5	Credits	23
5.1	Development	23
5.2	Supervision	23
5.3	Testing	23
5.4	External Code	23
5.5	Funding	23
5.6	Citation	24
6	Indices and tables	25
	Python Module Index	27

FISSA (Fast Image Signal Separation Analysis) is a Python library for decontaminating somatic signals from two-photon calcium imaging data. It can read images in tiff format and ROIs in zips as exported by ImageJ; or operate with numpy arrays directly, which can be produced by importing files stored in other formats.

For details of the algorithm, please see our [companion paper](#) published in Scientific Reports.

FISSA is compatible with both Python 2.7 and Python 3.5+. Using Python 3 is strongly encouraged, as Python 2 will no longer be [maintained starting January 2020](#).

FISSA has been tested on Ubuntu 17.04 and on Windows Windows 10 with the [Anaconda](#) distribution.

Documentation, including the full API, is available online at <https://fissa.readthedocs.io>.

If you encounter a specific problem please [open a new issue](#). For general discussion and help with installation or setup, please see the [Gitter chat](#).

1.1 Usage

A general tutorial on the use of FISSA can be found at: <https://rochefort-lab.github.io/fissa/examples/Basic%20usage.html>

An example workflow with another Python toolbox (SIMA): <https://rochefort-lab.github.io/fissa/examples/SIMA%20example.html>

An example workflow importing data exported from a MATLAB toolbox (cNMF): <https://rochefort-lab.github.io/fissa/examples/cNMF%20example.html>

These notebooks can also be run on your own machine. To do so, you will need to [download a copy of the repository](#), unzip it and browse to the [examples](#) directory. Then, start up a jupyter notebook server to run our notebooks. If you're new to jupyter notebooks, an approachable tutorial can be found at <https://www.datacamp.com/community/tutorials/tutorial-jupyter-notebook>.

1.2 Installation

1.2.1 Installation on Windows

Basic prerequisites

Download and install, in the following order:

- (for Python 2.7 only) Microsoft Visual C++ Compiler for Python 2.7: <https://www.microsoft.com/en-us/download/details.aspx?id=44266>
- Python 2.7 or 3.5+ (recommended) Anaconda as the Python environment, available from <https://www.anaconda.com/download/>.

Installing FISSA

Open `Anaconda Prompt.exe`, which can be found through the Windows start menu or search, and type or copy-paste (by right clicking) the following:

```
conda install -c conda-forge shapely tifffile
```

Then, install FISSA by running the command

```
pip install fissa
```

To test if FISSA has been installed, enter the command

```
python
```

to go into the Python environment. Then type

```
import fissa
```

If no errors show up, FISSA is now installed. You can leave Python by typing `exit()`.

If you want to use the interactive plotting from the notebooks, you should also install the HoloViews plotting toolbox, as follows

```
conda install -c ioam holoviews
```

See *usage* above for details on how to use FISSA.

1.2.2 Installation on Linux

Before installing FISSA, you will need to make sure you have all of its dependencies (and the dependencies of its dependencies) installed.

Here we will outline how to do all of these steps, assuming you already have both Python and pip installed. It is highly likely that your Linux distribution ships with these.

Dependencies of dependencies

- `scipy` requires a Fortran compiler and BLAS/LAPACK/ATLAS.
- `shapely` requires GEOS.

- `Pillow` $\geq 3.0.0$ effectively requires a JPEG library.

These packages can be installed on *Debian/Ubuntu* with the following shell commands.

```
sudo apt-get update
sudo apt-get install gfortran libopenblas-dev liblapack-dev libatlas-dev libatlas-
↳base-dev
sudo apt-get install libgeos-dev
sudo apt-get install libjpeg-dev
```

Installing FISSA

For normal usage of FISSA, you can install the latest release version on PyPI using pip:

```
pip install fissa
```

To also install fissa along with the dependencies required to run our sample notebooks (which include plots rendered with holoviews) you should run the following command:

```
pip install fissa['plotting']
```

Afterwards, you can test to see if FISSA is install by running the command

```
python
```

to start an interactive python session. Then run

```
import fissa
```

at the python command prompt.

If no errors show up, FISSA is now installed. You can leave the interactive python session with the `exit ()` command, or CTRL+D.

1.3 Folder Structure

1.3.1 continuous_integration/

Contains files necessary for deploying tests on continuous integration servers. Users should ignore this directory.

1.3.2 examples/

Contains example code. You can load the notebooks as `.ipynb` directly in GitHub, or on your system if you know how to use jupyter notebooks.

1.3.3 examples/exampleData/

Contains example data. It a zipfile with region of interests from ImageJ. It also contains three tiff stacks, which have been downsampled and cropped from full data from the Rochefort lab.

1.3.4 fissa/

Contains the toolbox.

1.3.5 fissa/tests/

Contains tests for the toolbox, which are run to ensure it will work as expected.

1.4 Citing FISSA

If you use FISSA for your research, please cite the following paper in any resulting publications:

S. W. Keemink, S. C. Lowe, J. M. P. Pakan, E. Dylida, M. C. W. van Rossum, and N. L. Rochefort. FISSA: A neuropil decontamination toolbox for calcium imaging signals, *Scientific Reports*, **8**(1):3493, 2018. DOI:[10.1038/s41598-018-21640-2](https://doi.org/10.1038/s41598-018-21640-2).

For your convenience, the FISSA package ships with a copy of this citation in bibtex format, available at citation.bib.

1.5 License

Unless otherwise stated in individual files, all code is Copyright (c) 2015, Sander Keemink, Scott Lowe, and Nathalie Rochefort. All rights reserved.

This program is free software; you can redistribute it and/or modify it under the terms of the GNU General Public License as published by the Free Software Foundation; either version 2 of the License, or (at your option) any later version.

This program is distributed in the hope that it will be useful, but WITHOUT ANY WARRANTY; without even the implied warranty of MERCHANTABILITY or FITNESS FOR A PARTICULAR PURPOSE. See the GNU General Public License for more details.

You should have received a copy of the GNU General Public License along with this program. If not, see <http://www.gnu.org/licenses/>.

2.1 fissa package

2.1.1 Submodules

fissa.ROI module

The functions below were adapted from the sima package <http://www.losonczylab.org/sima> version 1.3.0.

License

This file is Copyright (C) 2014 The Trustees of Columbia University in the City of New York.

This program is free software; you can redistribute it and/or modify it under the terms of the GNU General Public License as published by the Free Software Foundation; either version 2 of the License, or (at your option) any later version.

This program is distributed in the hope that it will be useful, but WITHOUT ANY WARRANTY; without even the implied warranty of MERCHANTABILITY or FITNESS FOR A PARTICULAR PURPOSE. See the GNU General Public License for more details.

You should have received a copy of the GNU General Public License along with this program. If not, see <http://www.gnu.org/licenses/>.

`fissa.ROI.poly2mask` (*polygons*, *im_size*)
Converts polygons to a sparse binary mask.

```
>>> from fissa.ROI import poly2mask
>>> poly1 = [[0,0], [0,1], [1,1], [1,0]]
>>> poly2 = [[0,1], [0,2], [2,2], [2,1]]
>>> mask = poly2mask([poly1, poly2], (3, 3))
>>> mask[0].todense()
matrix([[ True, False, False],
```

(continues on next page)

(continued from previous page)

```
[ True,  True, False],
 [False, False, False]], dtype=bool)
```

Parameters

- **polygons** (*sequence of coordinates or sequence of Polygons*) – A sequence of polygons where each is either a sequence of (x,y) or (x,y,z) coordinate pairs, an Nx2 or Nx3 numpy array, or a Polygon object.
- **im_size** (*tuple*) – Final size of the resulting mask

Returns **masks** – A list of sparse binary masks of the points contained within the polygons, one mask per plane. Each mask is in linked list sparse matrix format.

Return type list of sparse matrices

fissa.core module

Main user interface for FISSA.

Authors:

- Sander W Keemink (swkeemink@scimail.eu)
- Scott C Lowe

```
class fissa.core.Experiment (images, rois, folder, nRegions=4, expansion=1,
                             alpha=0.1, ncores_preparation=None, ncores_separation=None,
                             method='nmf', lowmemory_mode=False, datahandler_
                             dler_custom=None)
```

Bases: `object`

Does all the steps for FISSA.

```
__init__ (images, rois, folder, nRegions=4, expansion=1, alpha=0.1, ncores_preparation=None,
          ncores_separation=None, method='nmf', lowmemory_mode=False, datahandler_
          dler_custom=None)
```

Initialisation. Set the parameters for your Fissa instance.

Parameters

- **images** (*str or list*) – The raw recording data. Should be one of:
 - the path to a directory containing TIFF files (string),
 - an explicit list of TIFF files (list of strings),
 - a list of array_like data already loaded into memory, each shaped (*frames, y-coords, x-coords*).

Note that each TIFF/array is considered a single trial.

- **rois** (*str or list*) – The roi definitions. Should be one of:
 - the path to a directory containing ImageJ ZIP files (string),
 - the path of a single ImageJ ZIP file (string),
 - a list of ImageJ ZIP files (list of strings),
 - a list of arrays, each encoding a ROI polygons,
 - a list of lists of binary arrays, each representing a ROI mask.

This can either be a single roiset for all trials, or a different roiset for each trial.

- **folder** (*str*) – Output path to a directory in which the extracted data will be stored.
- **nRegions** (*int*, *optional*) – Number of neuropil regions to draw. Use a higher number for densely labelled tissue. Default is 4.
- **expansion** (*float*, *optional*) – Expansion factor for the neuropil region, relative to the ROI area. Default is 1. The total neuropil area will be $nRegions * expansion * area(ROI)$.
- **alpha** (*float*, *optional*) – Sparsity regularization weight for NMF algorithm. Set to zero to remove regularization. Default is 0.1. (Not used for ICA method.)
- **ncores_preparation** (*int*, *optional* (default: *None*)) – Sets the number of subprocesses to be used during the data preparation steps (ROI and subregions definitions, data extraction from tifs, etc.). If set to *None* (default), there will be as many subprocesses as there are threads or cores on the machine. Note that this behaviour can, especially for the data preparation step, be very memory-intensive.
- **ncores_separation** (*int*, *optional* (default: *None*)) – Same as *ncores_preparation*, but for the separation step. Note that this step requires less memory per subprocess, and hence can often be set higher than *ncores_preparation*.
- **method** ({'nmf', 'ica'}, *optional*) – Which blind source-separation method to use. Either 'nmf' for non-negative matrix factorization, or 'ica' for independent component analysis. Default (recommended) is 'nmf'.
- **lowmemory_mode** (*bool*, *optional*) – If *True*, FISSA will load TIFF files into memory frame-by-frame instead of holding the entire TIFF in memory at once. This option reduces the memory load, and may be necessary for very large inputs. Default is *False*.
- **datahandler_custom** (*object*, *optional*) – A custom datahandler for handling ROIs and calcium data can be given here. See datahandler.py (the default handler) for an example.

calc_deltaf (*freq*, *use_raw_f0=True*, *across_trials=True*)

Calculate deltaf/f0 for raw and result traces.

The results can be accessed as `self.deltaf_raw` and `self.deltaf_result`. `self.deltaf_raw` is only the ROI trace instead of the traces across all subregions.

Parameters

- **freq** (*float*) – Imaging frequency, in Hz.
- **use_raw_f0** (*bool*, *optional*) – If *True* (default), use an f0 estimate from the raw ROI trace for both raw and result traces. If *False*, use individual f0 estimates for each of the traces.
- **across_trials** (*bool*, *optional*) – If *True*, we estimate a single baseline f0 value across all trials. If *False*, each trial will have their own baseline f0, and df/f0 value will be relative to the trial-specific f0. Default is *True*.

save_to_matlab ()

Save the results to a matlab file.

Can be found in `folder/matlab.mat`.

This will give you a filename.mat file which if loaded in Matlab gives the following structs: ROIs, result, raw.

If $df/f0$ was calculated, these will also be stored as *df_result* and *df_raw*, which will have the same format as *result* and *raw*.

These can be interfaced with as follows, for cell 0, trial 0:

- *ROIs.cell0.trial0{1}* polygon for the ROI
- *ROIs.cell0.trial0{2}* polygon for first neuropil region
- *result.cell0.trial0(1,:)* final extracted cell signal
- *result.cell0.trial0(2,:)* contaminating signal
- *raw.cell0.trial0(1,:)* raw measured cell signal
- *raw.cell0.trial0(2,:)* raw signal from first neuropil region

separate (*redo_prep=False, redo_sep=False*)
 Separate all the trials with FISSA algorithm.

After running *separate*, data can be found as follows:

self.sep Raw separation output, without being matched. Signal *i* for a specific cell and trial can be found as *self.sep[cell][trial][i,:]*.

self.result Final output, in order of presence in cell ROI. Signal *i* for a specific cell and trial can be found at *self.result[cell][trial][i, :]*. Note that the ordering is such that *i = 0* is the signal most strongly present in the ROI, and subsequent entries are in diminishing order.

self.mixmat The mixing matrix (how to go between *self.separated* and *self.raw* from the *separation_prep()* function).

self.info Information about separation routine, iterations needed, etc.

Parameters

- **redo_prep** (*bool, optional*) – Whether to redo the preparation. Default is *False*. Note that if this is true, we set *redo_sep = True* as well.
- **redo_sep** (*bool, optional*) – Whether to redo the separation. Default is *False*. Note that this parameter is ignored if *redo_prep* is set to *True*.

separation_prep (*redo=False*)
 Prepare and extract the data to be separated.

For each trial, performs the following steps:

- Load in data as arrays
- Load in ROIs as masks
- Grow and separate ROIs to define neuropil regions
- Using neuropil and original ROI regions, extract traces from data

After running this you can access the raw data (i.e. pre-separation) as *self.raw* and *self.rois*. *self.raw* is a list of arrays. *self.raw[cell][trial]* gives you the traces of a specific cell and trial, across cell and neuropil regions. *self.roi_polys* is a list of lists of arrays. *self.roi_polys[cell][trial][region][0]* gives you the polygon for the region for a specific cell, trial and region. *region=0* is the cell, and *region>0* gives the different neuropil regions. For separateable masks, it is possible multiple outlines are found, which can be accessed as *self.roi_polys[cell][trial][region][i]*, where *i* is the outline index.

Parameters redo (*bool, optional*) – If *False*, we load previously prepared data when possible. If *True*, we re-run the preparation, even if it has previously been run. Default is *False*.

`fissa.core.extract_func(inputs)`

Extract data using multiprocessing.

Parameters `inputs` (*list*) – list of inputs

0. image array
1. the rois
2. number of neuropil regions
3. how much larger neuropil region should be then central ROI

Returns

- *dict* – Data across cells.
- *dict* – Polygons for each ROI.

`fissa.core.separate_func(inputs)`

Extraction function for multiprocessing.

Parameters `inputs` (*list*) – list of inputs

0. Array with signals to separate
1. Alpha input to `npil.separate`
2. Method
3. Current ROI number

Returns

- *numpy.ndarray* – The raw separated traces.
- *numpy.ndarray* – The separated traces matched to the primary signal.
- *numpy.ndarray* – Mixing matrix.
- *dict* – Metadata for the convergence result.

fissa.datahandler module

FISSA functions to handle image and roi objects and return the right format.

If a custom version of this file is used (which can be defined at the declaration of the core FISSA Experiment class), it should have the same functions as here, with the same inputs and outputs.

Authors:

- Sander W Keemink <swkeemink@scimail.eu>
- Scott C Lowe <scott.code.lowe@gmail.com>

`fissa.datahandler.extracttraces(data, masks)`

Extracts a temporal trace for each spatial mask.

Parameters

- **data** (*array_like*) – Data array as made by `image2array`. Should be shaped (*frames*, *y*, *x*).
- **masks** (*list of array_like*) – List of binary arrays.

Returns Trace for each mask. Shaped (*len(masks)*, *n_frames*).

Return type `numpy.ndarray`

`fissa.datahandler.getmean(data)`

Determine the mean image across all frames.

Parameters `data` (*array_like*) – Data array as made by `image2array`. Should be shaped (*frames*, *y*, *x*).

Returns *y* by *x* array for the mean values

Return type `numpy.ndarray`

`fissa.datahandler.image2array(image)`

Loads a TIFF image from disk.

Parameters `image` (*str* or *array_like*) – Either a path to a TIFF file, or *array_like* data.

Returns A 3D array containing the data, with dimensions corresponding to (*frames*, *y_coordinate*, *x_coordinate*).

Return type `numpy.ndarray`

`fissa.datahandler.rois2masks(rois, data)`

Take the object *rois* and returns it as a list of binary masks.

Parameters

- **rois** (*string* or *list of array_like*) – Either a string with imagej roi zip location, list of arrays encoding polygons, or list of binary arrays representing masks
- **data** (*array*) – Data array as made by `image2array`. Must be shaped (*frames*, *y*, *x*).

Returns List of binary arrays (i.e. masks)

Return type `list`

fissa.datahandler_framebyframe module

FISSA functions to handle image and roi objects and return the right format.

If a custom version of this file is used (which can be defined at the declaration of the core FISSA Experiment class), it should have the same functions as here, with the same inputs and outputs.

Authors:

- Sander W Keemink <swkeemink@scimail.eu>
- Scott C Lowe <scott.code.lowe@gmail.com>

`fissa.datahandler_framebyframe.extracttraces(data, masks)`

Get the traces for each mask in masks from data.

Parameters

- **data** (*PIL.Image*) – An open *PIL.Image* handle to a multi-frame TIFF image.
- **masks** (*list of array_like*) – List of binary arrays.

Returns Trace for each mask. Shaped (*len(masks)*, *n_frames*).

Return type `numpy.ndarray`

`fissa.datahandler_framebyframe.getmean(data)`

Determine the mean image across all frames.

Parameters `data` (*PIL.Image*) – An open *PIL.Image* handle to a multi-frame TIFF image.

Returns *y*-by-*x* array for the mean values.

Return type `numpy.ndarray`

`fissa.datahandler_framebyframe.image2array(image)`

Open a given image file as a `PIL.Image` instance.

Parameters `image` (*str or file*) – A filename (string) of a TIFF image file, a `pathlib.Path` object, or a file object.

Returns Handle from which frames can be loaded.

Return type `PIL.Image`

`fissa.datahandler_framebyframe.rois2masks(rois, data)`

Take the object 'rois' and returns it as a list of binary masks.

Parameters

- **rois** (*str or list of array_like*) – Either a string with imagej roi zip location, list of arrays encoding polygons, or list of binary arrays representing masks
- **data** (*PIL.Image*) – An open `PIL.Image` handle to a multi-frame TIFF image.

Returns List of binary arrays (i.e. masks).

Return type `list`

fissa.deltaf module

Functions for computing correcting fluorescence signals for changes in baseline activity.

Authors:

- Scott C Lowe

`fissa.deltaf.findBaselineF0(rawF, fs, axis=0, keepdims=False)`

Find the baseline for a fluorescence imaging trace line.

The baseline, F0, is the 5th-percentile of the 1Hz lowpass filtered signal.

Parameters

- **rawF** (*array_like*) – Raw fluorescence signal.
- **fs** (*float*) – Sampling frequency of rawF, in Hz.
- **axis** (*int, optional*) – Dimension which contains the time series. Default is 0.
- **keepdims** (*bool, optional*) – Whether to preserve the dimensionality of the input. Default is *False*.

Returns `baselineF0` – The baseline fluorescence of each recording, as an array.

Return type `numpy.ndarray`

Note: In typical usage, the input rawF is expected to be sized (*numROI, numTimePoints, numRecs*) and the output will then be sized (*numROI, 1, numRecs*) if *keepdims* is *True*.

fissa.neuropil module

Functions for removal of neuropil from calcium signals.

Authors:

- Sander W Keemink (swkeemink@scimail.eu)
- Scott C Lowe

Created: 2015-05-15

`fissa.neuropil.LowPassFilter(F, fs=40, nfilt=40, fw_base=10, axis=0)`

Low pass filters a fluorescence imaging trace line.

Parameters

- **F** (*array_like*) – Fluorescence signal.
- **fs** (*float, optional*) – Sampling frequency of F, in Hz. Default is 40.
- **nfilt** (*int, optional*) – Number of taps to use in FIR filter. Default is 40.
- **fw_base** (*float, optional*) – Cut-off frequency for lowpass filter, in Hz. Default is 10.
- **axis** (*int, optional*) – Along which axis to apply low pass filtering. Default is 0.

Returns Low pass filtered signal with the same shape as *F*.

Return type `numpy.ndarray`

`fissa.neuropil.separate(S, sep_method='nmf', n=None, maxiter=10000, tol=0.0001, random_state=892, maxtries=10, W0=None, H0=None, alpha=0.1)`

For the signals in *S*, finds the independent signals underlying it, using ica or nmf.

Parameters

- **S** (*array_like*) – 2-d array containing mixed input signals. Each column of *S* should be a different signal, and each row an observation of the signals. For $S[i,j]$, j = each signal, i = signal content. The first column, $j = 0$, is considered the primary signal and the one for which we will try to extract a decontaminated equivalent.
- **sep_method** (`{ 'ica', 'nmf' }`) – Which source separation method to use, either ICA or NMF.
 - `'ica'`: Independent Component Analysis
 - `'nmf'`: Non-negative Matrix Factorization
- **n** (*int, optional*) – How many components to estimate. If *None* (default), use PCA to estimate how many components would explain at least 99% of the variance and adopt this value for *n*.
- **maxiter** (*int, optional*) – Number of maximally allowed iterations. Default is 500.
- **tol** (*float, optional*) – Error tolerance for termination. Default is 1e-5.
- **random_state** (*int, optional*) – Initial state for the random number generator. Set to *None* to use the `numpy.random` default. Default seed is 892.
- **maxtries** (*int, optional*) – Maximum number of tries before algorithm should terminate. Default is 10.
- **W0** (*array_like, optional*) – Optional starting condition for *W* in NMF algorithm. (Ignored when using the ICA method.)
- **H0** (*array_like, optional*) – Optional starting condition for *H* in NMF algorithm. (Ignored when using the ICA method.)
- **alpha** (*float, optional*) – Sparsity regularization weight for NMF algorithm. Set to zero to remove regularization. Default is 0.1. (Ignored when using the ICA method.)

Returns

- `numpy.ndarray` – The raw separated traces.
- `numpy.ndarray` – The separated traces matched to the primary signal, in order of matching quality (see Methods below).
- `numpy.ndarray` – Mixing matrix.
- `dict` – Metadata for the convergence result, with keys:
 - `'random_state'`: seed for ICA initiation
 - `'iterations'`: number of iterations needed for convergence
 - `'max_iterations'`: maximum number of iterations allowed
 - `'converged'`: whether the algorithm converged or not (bool)

Notes

Concept by Scott Lowe and Sander Keemink. Normalize the columns in estimated mixing matrix A so that $\text{sum}(\text{column})=1$. This results in a relative score of how strongly each separated signal is represented in each ROI signal.

fissa.readimagejrois module

Based on code originally written by Luis Pedro Coelho <luis@luispedro.org>, 2012, available at <https://gist.github.com/luispedro/3437255>, distributed under the MIT License.

Modified

- 2014 by Jeffrey Zaremba (@jzaremba), <https://github.com/losonczylab/sima>
- 2015 by Scott Lowe (@scottclowe) and Sander Keemink (@swkeemink).

`fissa.readimagejrois.read_imagej_roi_zip(filename)`

Reads an ImageJ ROI zip set and parses each ROI individually

Parameters `filename` (*str*) – Path to the ImageJ ROIs zip file

Returns `roi_list` – List of the parsed ImageJ ROIs

Return type `list`

`fissa.readimagejrois.read_roi(roi_obj)`

Parses an individual ImageJ ROI

This is based on the Java implementation: <http://rsbweb.nih.gov/ij/developer/source/ij/io/RoiDecoder.java.html>

Parameters `roi_obj` (*file object*) – File object containing a single ImageJ ROI

Returns Returns a parsed ROI object, a dictionary with either a `'polygons'` or a `'mask'` field.

Return type `dict`

Raises

- `IOError` – If there is an error reading the roi file object
- `ValueError` – If unable to parse ROI

fissa.roitools module

Functions used for ROI manipulation.

Authors:

- Sander W Keemink <swkeemink@scimail.eu>

`fissa.roitools.find_roi_edge(mask)`

Finds the outline of a mask, using the `find_contour` function from `skimage.measure`.

Parameters `mask` (*array_like*) – the mask, a binary array

Returns `outline` – Array with coordinates of pixels in the outline of the mask

Return type list of (n,2)-ndarrays

`fissa.roitools.get_mask_com(mask)`

Get the center of mass for a boolean mask.

Parameters `mask` (*array_like*) – A two-dimensional boolean-mask.

Returns

- *float* – Center of mass along first dimension.
- *float* – Center of mass along second dimension.

`fissa.roitools.get_npil_mask(mask, totalexpansion=4)`

Given the masks for a ROI, find the surrounding neuropil.

Parameters

- **mask** (*array_like*) – The reference ROI mask to expand the neuropil from. The array should contain only boolean values.
- **expansion** (*float, optional*) – How much larger to make the neuropil total area than mask area.

Returns A boolean `numpy.ndarray` mask, where the region surrounding the input is now `True` and the region of the input mask is `False`.

Return type `numpy.ndarray`

Notes

Our implementation is as follows:

- On even iterations (where indexing begins at zero), expand the mask in each of the 4 cardinal directions.
- On odd numbered iterations, expand the mask in each of the 4 diagonal directions.

This procedure generates a neuropil whose shape is similar to the shape of the input ROI mask.

Note: For fixed number of *iterations*, squarer input masks will have larger output neuropil masks.

`fissa.roitools.getmasks(rois, shpe)`

Get the masks for the specified rois.

Parameters

- **rois** (*list*) – list of roi coordinates. Each roi coordinate should be a 2d-array or equivalent list. I.e.: `roi = [[0,0], [0,1], [1,1], [1,0]]` or `roi = np.array([[0,0], [0,1], [1,1], [1,0]])` I.e. a n by 2 array, where n is the number of coordinates. If a 2 by n array is given, this will be transposed.
- **shape** (*array/list*) – shape of underlying image [width,height]

Returns List of masks for each roi in the rois list

Return type list of numpy.ndarray

`fissa.roitools.getmasks_npil (cellMask, nNpil=4, expansion=1)`

Generate neuropil masks using the `get_npil_mask` function.

Parameters

- **cellMask** (*array_like*) – the cell mask (boolean 2d arrays)
- **nNpil** (*int*) – number of neuropil subregions
- **expansion** (*float*) – How much larger to make neuropil subregion area than in *cellMask*

Returns Returns a list with soma + neuropil masks (boolean 2d arrays)

Return type list

`fissa.roitools.readrois (roiset)`

read the imagej rois in the zipfile roiset, and make sure that the third dimension (i.e. frame number) is always zero.

Parameters **roiset** (*str*) – folder to a zip file with rois

Returns Returns the rois as polygons

Return type list

`fissa.roitools.shift_2d_array (a, shift=1, axis=None)`

Shifts an entire array in the direction of axis by the amount shift, without refilling the array.

Parameters

- **a** (*array_like*) – Input array.
- **shift** (*int, optional*) – How much to shift array by. Default is 1.
- **axis** (*int, optional*) – The axis along which elements are shifted. By default, the array is flattened before shifting, after which the original shape is restored.

Returns Array with the same shape as a, but shifted appropriately.

Return type numpy.ndarray

`fissa.roitools.split_npil (mask, centre, num_slices, adaptive_num=False)`

Splits a mask into a number of approximately equal slices by area around the center of the mask.

Parameters

- **mask** (*array_like*) – Mask as a 2d boolean array.
- **centre** (*tuple*) – The center co-ordinates around which the mask will be split.
- **num_slices** (*int*) – The number of slices into which the mask will be divided.
- **adaptive_num** (*bool, optional*) – If True, the *num_slices* input is treated as the number of slices to use if the ROI is surrounded by valid pixels, and automatically reduces the number of slices if it is on the boundary of the sampled region.

Returns A list with *num_slices* many masks, each of which is a 2d boolean numpy array.

Return type `list`

3.1 Reporting Issues

If you encounter a problem when implementing or using FISSA, we want to hear about it!

3.1.1 Gitter

To get help resolving implementation difficulties, or similar one-off problems, please ask for help on our [gitter channel](#).

3.1.2 Reporting Bugs and Issues

If you experience a bug, please report it by opening a [new issue](#). When reporting issues, please include code that reproduces the issue and whenever possible, an image that demonstrates the issue. The best reproductions are self-contained scripts with minimal dependencies.

Make sure you mention the following things:

- What did you do?
- What did you expect to happen?
- What actually happened?
- What versions of FISSA and Python are you using, and on which operating system?

3.1.3 Feature requests

If you have a new feature or enhancement to an existing feature you would like to see implemented, please check the list of [existing issues](#) and if you can't find it make a [new issue](#) to request it. If you do find it in the list, you can post a comment saying +1 (or :+1: if you are a fan of emoticons) to indicate your support for this feature.

3.2 Documentation

We are glad to accept any sort of documentation: function docstrings, tutorials, Jupyter notebooks demonstrating implementation details, etc.

reStructuredText documents and notebooks live in the source code repository under the `doc/` directory.

3.2.1 Docstrings

Documentation for classes and functions should follow the [format prescribed for numpy](#).

A complete example of this is available [here](#).

3.3 How to contribute

The preferred way to contribute to FISSA is to fork the [main repository](#) on GitHub.

1. Fork the [project repository](#). Click on the ‘Fork’ button near the top of the page. This creates a copy of the code under your account on the GitHub server.
2. Clone this copy to your local disk:

```
$ git clone git@github.com:YourUserName/fissa.git
$ cd fissa
```

3. Create a branch to hold and track your changes

```
$ git checkout -b my-feature
```

and start making changes.

4. Work on this copy on your computer using Git to do the version control. When you’re done editing, do:

```
$ git add modified_files
$ git commit
```

to record your changes in Git, writing a commit message following the [specifications below](#), then push them to GitHub with:

```
$ git push -u origin my-feature
```

5. Finally, go to the web page of your fork of the FISSA repo, and click ‘Pull request’ to issue a [pull request](#) and send your changes to the maintainers for review. This will also send a notification email to the committers.

If any of the above seems like magic to you, then look up the [Git documentation](#) on the web.

It is recommended to check that your contribution complies with the following rules before submitting a pull request.

- All public functions and methods should have informative docstrings, with sample usage included in the doctest format when appropriate.
- All unit tests pass. Check with (from the top level source folder):

```
$ pip install pytest
$ py.test
```

- Code with good unit test coverage (at least 90%, ideally 100%). Check with

```
$ pip install pytest pytest-cov
$ py.test --cov=fissa --cov-config .coveragerc
```

and look at the value in the ‘Cover’ column for any files you have added or amended.

If the coverage value is too low, you can inspect which lines are or are not being tested by generating a html report. After opening this, you can navigate to the appropriate module and see lines which were not covered, or were only partially covered. If necessary, you can do this as follows:

```
$ py.test --cov=fissa --cov-config .coveragerc --cov-report html --cov-report-
↳term-missing
$ sensible-browser ./htmlcov/index.html
```

- No **pyflakes** warnings. Check with:

```
$ pip install pyflakes
$ pyflakes path/to/module.py
```

- No **PEP8** warnings. Check with:

```
$ pip install pep8
$ pep8 path/to/module.py
```

AutoPEP8 can help you fix some of the easier PEP8 errors.

```
$ pip install autopep8
$ autopep8 -i -a -a path/to/module.py
```

Note that using the `-i` flag will modify your existing file in-place, so be sure to save any changes made in your editor beforehand.

These tests can be collectively performed in one line with:

```
$ pip install -r requirements-dev.txt
$ py.test --flake8 --cov=fissa --cov-config .coveragerc --cov-report html --cov-
↳report term
```

3.3.1 Commit messages

Commit messages should be clear, precise and stand-alone. Lines should not exceed 72 characters.

It is useful to indicate the nature of your commits with a commit flag, as described in the [numpy development guide](#).

You can use these flags at the start of your commit messages:

```
API: an (incompatible) API change
BLD: change related to building the package
BUG: bug fix
CI: change continuous integration build
DEP: deprecate something, or remove a deprecated object
DEV: development tool or utility
DOC: documentation; only change/add/remove docstrings, markdown or comments
ENH: enhancement; add a new feature without removing existing features
MAINT: maintenance commit (refactoring, typos, etc.); no functional change
REV: revert an earlier commit
RF: refactoring
STY: style fix (whitespace, PEP8)
```

(continues on next page)

(continued from previous page)

TST: addition or modification of tests REL: related to releases

3.4 Notes

This document was based on the contribution guidelines for [sklearn](#), [numpy](#) and [Pillow](#).

All notable changes to this project will be documented here.

The format is based on [Keep a Changelog](#), and this project adheres to [Semantic Versioning](#).

Categories for changes are: Added, Changed, Deprecated, Removed, Fixed, Security.

4.1 Version 0.6.1

Release date: 2019-03-11. Full commit changelog [on github](#).

4.1.1 Fixed

- Allow `deltaf.findBaselineF0` to run with fewer than 90 samples, by reducing the pad-length if necessary. (#62)
- Basic usage notebook wasn't supplying the correct `datahandler_custom` argument for the custom datahandler (it was using `datahandler` instead, which is incorrect; this was silently ignored previously but will now trigger an error). (#62)
- Use `ncores_preparation` for preparation step, not `ncores_separation`. (#59)
- Only use `ncores_separation` for separation step, not all cores. (#59)
- Allow both byte strings and unicode strings to be arguments of functions which require strings. Previously, byte strings were required on Python 2.7 and unicode strings on Python 3. (#60)

4.2 Version 0.6.0

Release date: 2019-02-26. Full commit changelog [on github](#).

4.2.1 Added

- Python 3 compatibility. (#33)
- Documentation generation, with Sphinx, Sphinx-autodoc, and Napoleon. (#38)

4.3 Version 0.5.3

Release date: 2019-02-18. Full commit changelog [on github](#).

4.3.1 Fixed

- Fix f0 detection with low sampling rates. (#27)

4.4 Version 0.5.2

Release date: 2018-03-07. Full commit changelog [on github](#).

4.4.1 Changed

- The default alpha value was changed from 0.2 to 0.1. (#20)

4.5 Version 0.5.1

Release date: 2018-01-10. Full commit changelog [on github](#).

4.5.1 Added

- Possibility to define custom datahandler script for other formats
- Added low memory mode option to load larger tiffs frame-by-frame (#14)
- Added option to use ICA instead of NMF (not recommended, but is a lot faster).
- Added the option for users to define a custom data and ROI loading script. (#13)

4.5.2 Fixed

- Fixed custom datahandler usage. (#14)
- Documentation fixes. (#12)

4.6 Version 0.5.0

Release date: 2017-10-05

Initial release

5.1 Development

- Sander Keemink : Conception, Overall development
- Scott Lowe : Overall development

5.2 Supervision

- Nathalie Rochefort
- Mark van Rossum

5.3 Testing

- Janelle Pakan
- Evelyn Dylida

5.4 External Code

- `ROI.py` was adapted from code contained in the [SIMA](#) package under the [GNU GPL v2 License](#).
- `NIH ImageJ ROI parsing` was adapted from code originally written by Luis Pedro Coelho under the [MIT license](#).

5.5 Funding

- SK: DTC & Eurospin

- Wellcome Trust (Sir Henry Dale Fellowship), EU funding (Career Integration Grant)
- JP: EU intra-European Fellowship
- Centre for Integrative Physiology, University of Edinburgh

5.6 Citation

S. W. Keemink, S. C. Lowe, J. M. P. Pagan, E. Dylida, M. C. W. van Rossum, and N. L. Rochefort. FISSA: A neuropil decontamination toolbox for calcium imaging signals, *Scientific Reports*, **8**(1):3493, 2018. DOI:[10.1038/s41598-018-21640-2](https://doi.org/10.1038/s41598-018-21640-2).

CHAPTER 6

Indices and tables

- `genindex`
- `modindex`
- `search`

f

- `fissa`, [5](#)
- `fissa.core`, [6](#)
- `fissa.datahandler`, [9](#)
- `fissa.datahandler_framebyframe`, [10](#)
- `fissa.deltaf`, [11](#)
- `fissa.neuropil`, [11](#)
- `fissa.readimagejrois`, [13](#)
- `fissa.ROI`, [5](#)
- `fissa.roitools`, [14](#)

Symbols

`__init__()` (fissa.core.Experiment method), 6

C

`calc_deltaf()` (fissa.core.Experiment method), 7

E

`Experiment` (class in fissa.core), 6

`extract_func()` (in module fissa.core), 8

`extracttraces()` (in module fissa.datahandler), 9

`extracttraces()` (in fissa.datahandler_framebyframe), 10

F

`find_roi_edge()` (in module fissa.roitools), 14

`findBaselineF0()` (in module fissa.deltaf), 11

`fissa` (module), 5

`fissa.core` (module), 6

`fissa.datahandler` (module), 9

`fissa.datahandler_framebyframe` (module), 10

`fissa.deltaf` (module), 11

`fissa.neuropil` (module), 11

`fissa.readimagejrois` (module), 13

`fissa.ROI` (module), 5

`fissa.roitools` (module), 14

G

`get_mask_com()` (in module fissa.roitools), 14

`get_npil_mask()` (in module fissa.roitools), 14

`getmasks()` (in module fissa.roitools), 14

`getmasks_npil()` (in module fissa.roitools), 15

`getmean()` (in module fissa.datahandler), 10

`getmean()` (in module fissa.datahandler_framebyframe), 10

I

`image2array()` (in module fissa.datahandler), 10

`image2array()` (in fissa.datahandler_framebyframe), 11

L

`lowPassFilter()` (in module fissa.neuropil), 12

P

`poly2mask()` (in module fissa.ROI), 5

R

`read_imagej_roi_zip()` (in module fissa.readimagejrois), 13

`read_roi()` (in module fissa.readimagejrois), 13

`readrois()` (in module fissa.roitools), 15

`rois2masks()` (in module fissa.datahandler), 10

`rois2masks()` (in fissa.datahandler_framebyframe), 11

S

`save_to_matlab()` (fissa.core.Experiment method), 7

`separate()` (fissa.core.Experiment method), 8

`separate()` (in module fissa.neuropil), 12

`separate_func()` (in module fissa.core), 9

`separation_prep()` (fissa.core.Experiment method), 8

`shift_2d_array()` (in module fissa.roitools), 15

`split_npil()` (in module fissa.roitools), 15